

List of PCR primers used in this study

<b>PCR amplification</b>	<b>Forward (5'-3')</b>	<b>Reverse (5'-3')</b>
5' homology arm	CAGGGTTATACACACAGAGGCCACAAGAAG	CCATGGCGCGGGGCGGGGTGACTTCG
3' homology arm	ATGAAAAGGAAGGTAGATTCACATATATAC	TTGATTACAGCCATGTTCTGATCTTG
5' external probe	TCTCTCCACCTGCCCAACTTGCATCA	CTAGTTCACATCAGGAGTCTCAAGACTGTG
3' external probe	GCTCTACCCTGGAGATTACATAGATG	GTAAGGTCTAAGGCATTGTCTTCAAAG
5' cDNA probe	ATACG TTCACGAGGGGTCCTGTTTG	CTGGTCAGTTGTCATCTGTACTGCC
3' cDNA probe	CAGCCACATGCAAAGTGCCTTTGTG	GTCATTGGCGTCTGTTTTCTTAGTGC
Real-time PCR for 3' mouse <i>Hltf</i> transcript	TGCGCAGCTTCTCGAGTGTTCTTA	ATCACCTCTTGCTTCTGGCCAAGT